



# Moveable Headwalls

Prepare for the unexpected with our Fast Track Moveable Axiom Headwall.

When responding to a crisis, you want flexibility and speed on your side. By simply rolling our moveable headwalls into place, these double-sided units provide immediate access to medical gas services for two patients.

We have streamlined the process for HSI's Fast Track Moveable Axiom Headwall, allowing for a more affordable price and faster turnaround time. With special financing, you put down a deposit on the units you might need and pay the balance only when your headwalls are needed onsite.

HSI's Fast Track Moveable Axiom Headwalls are designed to help you to prepare for the unexpected.



750 Garcia Avenue    Tel: 925.427.7800    ◇    info@HSIheadwalls.com  
Pittsburg, CA 94565    Fax: 925.427.0800    ◇    www.HSIheadwalls.com



Made in the USA

## Movable Axiom™ EZ Headwall

The **Moveable Axiom™ EZ Headwalls** are specifically designed for pandemic applications and offer a quick turn and submittal free, easy to connect option. The vertically configured headwalls provide Vacuum, Air and Oxygen with hose connections to increase available services to two (2) patients in emergency situations. All hose connections are located on the right hand side of the unit to allow for maximum available movement from the building wall. All power is live by plugging the unit into an electrical receptacle using the 6ft power cord.

- Moveable Units are available to lease for a period no less than 4 weeks and for as long as required by facility
- When units are ready to be returned to HSI, customers can re-use the crates provided, or HSI can build new crates and ship to site for an additional cost
- All Vacuum gas outlets will be required to be capped using HSI provided materials prior to shipment back to HSI
- All Medical Gas hoses are required to be disposed of by the facility or customer, (HSI will not take back any hoses)
- All leased units are treated with MicrobeCare™ Non leaching, instant kill, Anti-Microbial application



### Features

- Locking Casters, 2 Fixed and 2 Moveable, for easy maneuverability

NOTE:  
 CONNECTIONS BETWEEN UNIT AND IN WALL MEDICAL GAS OUTLETS ARE MADE VIA MEDICAL GAS HOSE

BEACON MADEAS MEDICAL GAS OUTLETS	
KEY STYLE:	[ ] PURITAN BENNETT (QC) [ ] CHEMETRON (QC) [ ] MEDAES (OHMEDA) (QC) [ ] DISS
PIPE SIZES "TYPE L"	REGULATORY AGENCY APPROVAL
OXY 1/2" NOM.	U.L. LISTED
AIR 1/2" NOM.	
VAC 3/4" NOM. ASTM B819	

QTY. REQ'D.: \_\_\_\_\_ UNIT

CEILING HEIGHT:  
 \_\_\_\_\_

MATERIAL:

FRAME ANODIZED ALUMINUM  
 FASCIA COLD ROLLED STEEL

FINISH:

FRAME DEEP ETCHED, CLEAR ANODIZED  
 FASCIA POWDER COATED OYSTER WHITE

**GENERAL NOTES**

BASIC FRAME IS CONSTRUCTED OF DEEP ETCHED, CLEAR ANODIZED EXTRUDED ALUMINUM.

UNIMOUNT™ ACCESSORY TRACKS ARE INTEGRAL ON BOTH SIDES OF THE UNIT.




FASCIA ARE CONSTRUCTED OF COLD ROLLED STEEL POWDER COATED OYSTER WHITE, AND ARE HELD IN PLACE WITH MACHINE SCREWS CONCEALED BY GRAY PLASTIC VINYL TRIM.

DEVICE COVER PLATES ARE CLEAR ANODIZED ALUMINUM. MEDICAL GAS COVER PLATES ARE PER MEDICAL GAS MANUFACTURER'S STANDARD.

MEDICAL GAS SERVICE CONNECTIONS AND ELECTRICAL/COMMUNICATION CONNECTIONS SHALL BE THROUGH THE TOP OF THE UNIT. PRE-MANIFOLDED MEDICAL GAS PIPING WILL TERMINATE AS INDICATED ON THE DRAWING.

**LEGEND**

ITEM	DESCRIPTION
O	OXYGEN OUTLET
A	AIR OUTLET
V	VACUUM OUTLET
N	RECEPTACLE, DUPLEX, NORMAL, WHITE



**HSI**  
 Hospital Systems, Inc.  
 750 Garcia Avenue, Pittsburg, CA, 94565 USA  
 TEL: +1 925 427 7800 - FAX: +1 925 427 0800  
 www.hsiheadwalls.com

THIS DOCUMENT WITH THE INFORMATION PRINTED HEREON ARE CONFIDENTIAL PROPRIETARY PROPERTY OF HOSPITAL SYSTEMS, INC. AND ARE FURNISHED SOLELY FOR USE IN OBTAINING APPROVAL OF THE CONSTRUCTION DETAILS AND MATERIALS. NO PART OF THIS DOCUMENT MAY BE REPRODUCED, USED OR DISTRIBUTED EITHER IN WHOLE OR IN PART WITHOUT AUTHORIZATION FROM HOSPITAL SYSTEMS, INC.

REVISION	
1	
2	
3	
4	
5	

TITLE: 2-SIDED MOVABLE AXIOM HEADWALL  
 PROJECT: HSI HEADWALLS  
 MOVEABLE AXIOM EZ, TRIAGE

Drawn:	SR
Date:	200327
Checked:	
Scale:	1/16":1"

Part Number:  
**1901845**



MICROBECARE™

**HSI**  
Hospital Systems, Inc.  
*The Headwall Experts*

MicrobeCare™ is a unique antimicrobial bonding technology that effectively controls bacteria, fungi, algae, and yeast on a wide variety of treated materials and substrates. It has been shown to kill > 99.99 percent of microbes and does not allow for microbe mutation. The antimicrobial active agent in MicrobeCare™ is registered with the US Environmental Protection Agency (EPA). The antimicrobial coating has been used safely and effectively and has recorded proven results in highly sensitive healthcare environments.

### How does MicrobeCare™ work?

- A covalent bond permanently affixes MicrobeCare™ to the surface
- Cross linking between adjacent molecules forms a matrix of swords
- The positively charged nitrogen atom attracts microbes to the surface and ultimately electrocutes the microorganisms through a positive—negative ion exchange
- 18 link carbon chain pierces microbes as they are pulled towards the nitrogen atom

### Non-Leaching

Most antimicrobials use leaching technologies which offer a “zone of inhibition” unfortunately, leaching allows the antimicrobial to migrate into the environment where no antimicrobial was applied. As the antimicrobial leaches further from the application site, it exists in a weakened state which will eventually allow for cohabitation of microorganisms.

### Patient Room 'High Touch' areas



HSI is honored to be the sole licensee of Parasol Medical to apply MicrobeCare™ to our Headwalls and Accessories.

Any questions please reach out.

750 Garcia Avenue    Tel: 925.427.7800    ◇    [info@HSIheadwalls.com](mailto:info@HSIheadwalls.com)  
Pittsburg, CA 94565    Fax: 925.427.0800    ◇    [www.HSIheadwalls.com](http://www.HSIheadwalls.com)